



Spare Parts Manual

1. MilkoScan™ FT3 — Spare Parts Manual

2. Introduction

2.1. Maintenance Procedures

2.1.1. 1. Pre-verification of Instrument before PM replacement

2.1.2. 2. Removal and Replacement

2.1.3. 3. Post-verification of Instrument after PM

2.1.4. 4. Instrument Check

2.1.5. 5. Software

2.1.6. 6. Instrument Performance Verification

3. Maintenance Procedures (w/Index)

3.1. 1. Pre-verification of Instrument before PM replacement

3.2. 2. Removal and Replacement

3.3. 3. Post-verification of Instrument after PM

3.4. 4. Instrument Check

3.5. 5. Software

3.6. 6. Instrument Performance Verification

4. Untitled

1. MilkoScan™ FT3 — Spare Parts Manual

Document No.: 60085771 / Rev. 4

Copyright: 2020 / All rights reserved

Publisher: FOSS Analytical A/S, Nils Foss Allé 1, DK-3400 Hillerød, Denmark

Tel: +45 7010 3370 | **Fax:** +45 7010 3371

All information is liable to change without prior notice. For latest information about documentation updates for your specific instrument, please contact your local FOSS representative.

Revision History

| Rev. | Date of Issue | Revised Material |
|------|---------------|---|
| 1 | 2020-04-09 | First issue. |
| 2 | 2020-08-19 | Major update of spare parts and drawings. |
| 3 | 2022-01-24 | Major update of spare parts and drawings. |
| 4 | 2022-06-22 | Major update of spare parts and drawings. |

Table of Contents

1. General

1.1 Explanation of Definitions

2. Waste Sensor

3. Pre-heater

4. Sample Pipette - Sample Intake Mod

5. [Sample Pipette Mod - Main Assy](#)
 6. [Sample Pipette Mod - Main Assy \(cont.\)](#)
 7. [System Pipette Mod](#)
 8. [IR-source](#)
 9. [Filterstick for IFU Manifold 24 µm](#)
 10. [Filterstick for IFU Manifold 34 µm](#)
 11. [Detector](#)
 12. [IF Unit Main](#)
 13. [IF Unit Main \(cont.\)](#)
 14. [BPV Mod Main](#)
 15. [M2 Pump Mod - Main Assy](#)
 16. [H-stage](#)
 17. [M1 Pump Mod - Main Assy](#)
 18. [CTR Mod - Main Assy](#)
 19. [Main Unit](#)
 20. [Main Unit \(cont.\)](#)
 21. [Main Unit \(cont.\)](#)
 22. [Main Unit - Tube Diagram](#)
 23. [Main Unit - Connection Diagram](#)
 24. [Basic Unit with Conductivity Module](#)
 25. [Index](#)
-

1. General

[Go to Part Index](#)

This document contains Spare Parts for MilkoScan™ FT3.

Warning: Please read Safety Precautions in the Service Manual prior to maintenance or exchange of spare parts. The system may only be repaired by personnel certified by FOSS.

1.1. Explanation of Definitions

Explanation of definitions used in this manual:

| Term | Definition |
|-------------|---------------------------------|
| Drawing No. | FOSS drawing number. |
| Page | Drawing page in FOSS database. |
| Part No. | Spare part number for ordering. |
| Prod. No. | FOSS production number. |

Note: When ordering new spare parts, please use the Part No. stated in this manual for the parts in question.

2. Waste Sensor

Drawing No.: 60093513 | **Page:** 100-1

| Pos. | Part No. | Description | Prod. No. |
|------|----------|---------------------------|-----------|
| 2 | 60085248 | BNC Male to double banana | |
| 4 | 60078292 | Tube ØØ008/016 Silicone | |

(See Fig. 1)





Fig. 1

3. Pre-heater

Drawing No.: 60086467 | **Page:** 110-1

| Pos. | Part No. | Description | Prod. No. |
|------|----------|-------------------------------|-----------|
| 4 | 60081085 | Heating Cartridge and NTC cpl | |

(See Fig. 2)



600647 Doc. Rev 01 Doc Type 70




Fig. 2


Fig. 2

4. Sample Pipette - Sample Intake Mod

Drawing No.: 60073632 | Page: 200-1

| Pos. | Part No. | Description | Prod. No. |
|------|----------|--|-----------|
| 3 | 60071292 | Pipette Light Guide - Sample Pipette Mod | |

(See Fig. 3)



60073632 Doc. Rev 10 DocType95

200

Fig. 3


Fig. 3

5. Sample Pipette Mod - Main Assy

Drawing No.: 60071249 | Page: 210-1

| Pos. | Part No. | Description | Prod. No. |
|------|----------|--|-----------|
| 2 | 60079612 | Sensor optical switch fork for pipette c | |
| 3 | 60064537 | Miniature Speaker 8 Ohm, wire L=180mm | |
| 5 | 60040345 | Lid dampers | |
| 15 | 00208538 | O-ring | |
| 18 | 60071294 | Shaft For Pipette Rotation - Sample Inta | |
| 21 | 60034902 | Vibrationdamper | |
| 31 | 60075728 | Spur Gear m=0.8, for pipette rot mech | |
| 32 | 60080240 | Pipette Motor Cpl. | |

(See Fig. 4)



60079612 Doc. Rev. 10 Dec Page 10

210-1

Fig. 4

Fig. 4

6. Sample Pipette Mod - Main Assy (cont.)

Drawing No.: 60071249 | Page: 210-2

| Pos. | Part No. | Description | Prod. No. |
|------|----------|--|-----------|
| 26 | 60093006 | Sample Pipette Assy - Sample Intake Mod | 60073632 |
| 30 | 60077361 | Waste Funnel - Sample Intake | |
| A | 60073182 | Sample holder platform with holes placed | |

Note: For installation this part is located in 60084298 FOSS Inst.Kit. MilkoScan FT3.

(See Fig. 5)




Fig. 5

Fig. 5

7. System Pipette Mod

| Pos. | Part No. | Description | Prod. No. |
|------|----------|---|-----------|
| 4 | 60073622 | Liquid Sensor Cpl - System Pipette Assy | |
| 6 | 00039388 | Spring | |
| 500 | 60077482 | Cable ERNI SMC 12Pin 300mm | |
| 501 | 60080569 | Cable ERNI PVC (awg26) 4pin 300mm | |
| 502 | 60088654 | Cable Blade 6.3 F/Ring M3 1x22 L100 blk | |

(See Fig. 6)





Fig. 6

Fig. 6

8. IR-source

| Pos. | Part No. | Description | Prod. No. |
|------|----------|-----------------------------------|-----------|
| 4 | 60073905 | IR source Irma | |
| 8 | 00209452 | O-ring | |
| 500 | 60071333 | Cable Erni SMC 12P F/F 12×30 L150 | |

(See Fig. 7)



60073908 Doc. Rev 10 DocType 95


Fig. 7

Fig. 7

9. Filterstick for IFU Manifold 24 µm

| Pos. | Part No. | Description | Prod. No. |
|------|----------|----------------------|-----------|
| 2 | 60015107 | O-ring Ø003.00xØ1.00 | |

(See Fig. 8)



60070544 Doc. Rev.05 Doc. type:70

520

Fig. 8


Fig. 8

10. Filterstick for IFU Manifold 34 µm

Drawing No.: 60070544 | Page: 530-1

| Pos. | Part No. | Description | Prod. No. |
|------|----------|----------------------|-----------|
| 2 | 60015107 | O-ring Ø003.00xØ1.00 | |

(See Fig. 9)



6007044 Doc. Rev. 01 Doc. Type 3D

530

Fig. 9


Fig. 9

11. Detector

Drawing No.: 60070995 | Page: 540-1

| Pos. | Part No. | Description | Prod. No. |
|------|----------|-----------------------------------|-----------|
| 1 | 60093007 | Detector cpl | 60076967 |
| 5 | 60070130 | Detector PCB Lid | |
| 500 | 60077478 | Cable Erni SMC 26P F/F 26×26 L160 | |

(See Fig. 10)



60071578 Doc. Rev. 02 DocType 05

Fig. 10


Fig. 10

12. IF Unit Main

Drawing No.: 60071578 | Page: 550-1

| Pos. | Part No. | Description | Prod. No. |
|----------|----------|-----------------------------------|-----------|
| 19 | 60070021 | Filterstick for IFU Manifold 34µm | 60070544 |
| 20 | 60070134 | Fingerscrew for Ø4 | |
| 26 | 60084455 | IFU manifold | |
| 31 | 60072743 | Bulkhead IFU Gasket - Main Assy | |
| 32 | 60093008 | IR-source assy | 60073908 |
| 33 | 60084840 | LEDA2 IFU | 60072810 |
| 33-4a | 01052263 | Pcb assy IFU Control PCB | |
| 33-14-21 | 01052305 | Pcb assy Thermobox Heater PCB | |
| 35 | 60092630 | Cuvette MilkoScan FT3 assy | 60084456 |

(See Fig. 11)



60070544 Doc. Rev 20 Doc. Type 95

Fig. 11


Fig. 11

13. IF Unit Main (cont.)

Drawing No.: 60071578 | Page: 550-2

| Pos. | Part No. | Description | Prod. No. |
|------|----------|--------------------------------|-----------|
| 28 | 60075011 | Damper Lord mount | |
| 33-6 | 60078531 | IFU Dryer Cable | |
| 500 | 60080871 | Sensor optical switch fork cpl | |

(See Fig. 12)



60071578 Doc. Rev.20 Doc.Type:PS

550-2

Fig. 12

14. BPV Mod Main

| Pos. | Part No. | Description | Prod. No. |
|------|----------|-----------------------------------|-----------|
| 4 | 60071155 | Valve solenoid 2/2 NC 24V Bürkert | |
| 6 | 60053831 | Membrane | |
| 11 | 60087522 | BPV Manifold Cpl | |
| A | 60093500 | Motor step Kit for BPV | |

(See Fig. 13)




Fig. 13

Fig. 13

15. M2 Pump Mod - Main Assy

| Pos. | Part No. | Description | Prod. No. |
|------|----------|------------------------------------|-----------|
| 15 | 60085631 | Inner Housing w/ball - H-stage Cpl | |
| A | 60093501 | Motor step Kit for H-stage | |

(See Fig. 15)




Fig. 15

Fig. 15

17. M1 Pump Mod - Main Assy


Drawing No.: 60071426 | Page: 810-1

| Pos. | Part No. | Description | Prod. No. |
|------|----------|--|-----------|
| 8 | 60091950 | H-stage | 60085620 |
| 10 | 60093009 | Pressure Sensor Ø11, 0-30 bar (abs) Cpl | 60077108 |
| 11 | 60079419 | Heating cartridge Ø6.5×50mm 24V/55W cpl | |
| 12 | 60090139 | Motor step 43×43mm, Screw L=2 in cpl | |
| 15 | 01052233 | Pcb assy Pipette Light PCB | |
| 18 | 60015423 | O-ring Ø002.50xØ1.00 | |
| 21 | 60050165 | O-ring Ø011.50xØ1.00 EPDM | |
| 22 | 60086575 | Valve solenoid 2/2 NC 24V Bürkert Ø0.8 H | |
| 24 | 60044049 | O-ring Ø004.00xØ2.00 Viton 75 sh | |
| 27 | 60083598 | Pressure Sensor Retainer | |
| 500 | 60078475 | Cable manifold light | |
| A | 60092440 | Repair kit for piston pump | |

(See Fig. 16)

| Pos. | Part No. | Description | Prod. No. |
|------|----------|---|-----------|
| 1 | 60078314 | Fan Axial 92×92×25 3414 cpl | |
| 7 | 60090830 | Power Supply | 60088699 |
| 23 | 01052407 | Pcb assy NAT ZYNQ PCB | |
| 26 | 01052347 | Pcb Control PCB | |
| 500 | 60078537 | Power cable 24V Control PCB | |
| 501 | 60078516 | Cable Blade 6.3 F/Ring M3 1×16 L130 y/g | |
| 502 | 60078517 | Cable Blade 6.3 F/F 1×20 L200 y/g | |
| 503 | 60078508 | Cable Blade 6.3 F/Ring M3 1×18 L130 Brn | |
| 504 | 60078519 | Cable Blade 6.3 F/F 1×18 L170 blu | |
| 505 | 60078530 | Cable Blade 6.3 F/F 1×18 L170 Brn | |
| 506 | 60078507 | Cable Blade 6.3 F/Ring M3 1×18 L130 Blu | |

(See Fig. 17)



60070150 Dec. Rev 05 Dec Page 05

Fig. 17

Fig. 17

19. Main Unit

Drawing No.: 60070150 | Page: 1000-1

| Pos. | Part No. | Description | Prod. No. |
|------|----------|---|-----------|
| 1-7 | 60084479 | Flange Head Bolt, M4×20, ISO 7380-2, A2 | |
| 1-9 | 00182782 | Washer | |
| 1-10 | 60078539 | Gas Spring Ø6 - Stabilus | |
| 1-12 | 60079472 | Damper Foot | |
| 1-19 | 60071790 | Front Hatch Cpl - Basic Cabinet Assy | 60045881 |
| 1-20 | 60073494 | Top Hatch Cpl - Basic Cabinet Assy | |
| 1-21 | 60078383 | Gas Spring Connector | |
| 1-22 | 60090207 | Flange Head Bolt, M4×8 ISO 7380-2 A2 | |
| 1-23 | 60076302 | Top lid seal short | |
| 7 | 60091949 | Display Mod - Main Assy | 60073326 |
| 8 | 60070022 | System Pipette Mod | 60073197 |
| 10 | 60086467 | Pre-heater assy | |
| 12 | 60078314 | Fan Axial 92×92×25 3414 cpl | |
| 48 | 60065855 | Holder for tube Ø4.0 12 pos. PP black | |
| 703 | 60075782 | Label Clean Liquid | |
| 704 | 60075783 | Label Zero Liquid | |

(See Fig. 18)




Fig. 18

Fig. 18

20. Main Unit (cont.)

Drawing No.: 60070150 | Page: 1000-2

| Pos. | Part No. | Description | Prod. No. |
|------|----------|----------------------------------|-----------|
| 13 | 60078510 | Fan Centrifugal 135×135×38 24VDC | |
| 15 | 60078940 | Damper, pneumatic | |
| 28 | 01052349 | Pcb assy Driver Connection PCB | 00584698 |
| 29 | 60048004 | Chassis ID | |

(See Fig. 19)




Fig. 19

Fig. 19

21. Main Unit (cont.)

Drawing No.: 60070150 | Page: 1000-3

| Pos. | Part No. | Description | Prod. No. |
|------|----------|------------------------------------|-----------|
| 500 | 60048020 | Cable, Chas-ID PCB L=90 #2×0,25 UL | 60023887 |

(See Fig. 20)




Fig. 20


Fig. 20

22. Main Unit - Tube Diagram

Drawing No.: 60070150 | Page: 1000-4

| Pos. | Part No. | Description | Prod. No. |
|------|----------|-------------------------------------|-----------|
| 26 | 60091427 | Tube Ø01.59/Ø03.18 PTFE L=0145 assy | |
| 35 | 60091424 | Tube Ø01.00/Ø01.59 PTFE L=075 assy | |
| 50 | 60086446 | Tube Ø01.59/Ø03.18 PTFE L=0300 assy | |
| 51 | 60088416 | Tube Ø01.00/Ø01.59 PTFE L=0130 assy | |
| 52 | 60074456 | Tube Ø01.00/Ø01.59 PTFE L=0090 assy | |
| 53 | 60074464 | Tube Ø01.59/Ø03.18 PTFE L=0725 assy | |
| 54 | 60075536 | Tube Ø01.59/Ø03.18 PTFE L=0415 assy | |
| 55 | 60074458 | Tube Ø01.59/Ø03.18 PTFE L=0255 assy | |
| 56 | 60070955 | Tube Ø01.59/Ø03.18 PTFE L=0130 assy | |

(See Fig. 21)



6007058 Doc. Rev 17 Doc type 05

1000-L

Fig. 21

Fig. 21

22.1 Tube Dimensions and Materials Reference

| ID | Dimensions and Material | Material # |
|----|---------------------------|------------|
| B | Ø1.59/3.18 × 130 PTFE | 60070955 |
| C | Ø1.00/1.59 × 105 PTFE | 60074456 |
| D | Ø1.59/3.18 × 300 PTFE | 60086446 |
| E | Ø1.59/3.18 × 415 PTFE | 60075536 |
| F | Ø1.00/1.59 × 130 PTFE | 60088416 |
| G | Ø1.59/3.18 × 725 PTFE | 60074464 |
| GC | Ø1.59/3.18 × 145 PTFE | 60091427 |
| H | Ø1.00/1.59 × 105 PTFE | 60074456 |
| I | Ø1.59/3.18 × 255 PTFE | 60074458 |
| J | Ø1.00/1.59 × 075 PTFE | 60091424 |
| K | Ø1.59/3.18 × 415 PTFE | 60075536 |
| M | Ø1.59/3.18 × 415 PTFE | 60075536 |
| Z | Ø8.0/Ø16.0 x 850 Silicone | 60078292 |

(See Fig. 24)




Fig. 24

Fig. 24

23. Main Unit - Connection Diagram

Drawing No.: 60070150 | Page: 1000-6

| Pos. | Part No. | Description | Prod. No. |
|------|----------|--|-----------|
| 501 | 60081489 | Cable Power Control/Driver L=0800 8 UL | |
| 502 | 60093452 | Ethernet cable 1:1 SF/FTP CAT.6 CMX | |
| 503 | 60081519 | Cable Minibridge 8P F/F 8×26 L620 | |

(See Fig. 23)




Fig. 23

Fig. 23

24. Basic Unit with Conductivity Module

Drawing No.: 60093239 | Page: 1100-1

| Pos. | Part No. | Description | Prod. No. |
|------|----------|-----------------------------|-----------|
| 2 | 60092766 | FT3 Conductivity Module bdl | 60092345 |

(See Fig. 25)




Fig. 25

Fig. 25

25. Index

| Part No. | Page(s) |
|----------|---------|
| 00039388 | 2, 4 |
| 00182782 | 36 |
| 00208538 | 8 |
| 00539056 | 12 |
| 01052233 | 32 |
| 01052263 | 22 |
| 01052305 | 22 |
| 01052347 | 34 |
| 01052349 | 38 |
| 01052407 | 34 |
| 60015107 | 16, 18 |
| 60015423 | 32 |
| 60034902 | 8 |
| 60040345 | 8 |
| 60044049 | 32 |
| 60048004 | 38 |
| 60048020 | 40 |
| 60050165 | 32 |
| 60053831 | 26 |
| 60064537 | 8 |
| 60065855 | 36 |
| 60070021 | 22 |
| 60070022 | 36 |
| 60070130 | 20 |
| 60070134 | 22 |

| Part No. | Page(s) |
|----------|---------|
| 60070955 | 42 |
| 60071155 | 26, 28 |
| 60071292 | 6 |
| 60071294 | 8 |
| 60071333 | 14 |
| 60071790 | 36 |
| 60072743 | 22 |
| 60073182 | 10 |
| 60073494 | 36 |
| 60073905 | 14 |
| 60074456 | 42 |
| 60074458 | 42 |
| 60074464 | 42 |
| 60075011 | 24 |
| 60075536 | 42 |
| 60075728 | 8 |
| 60075782 | 36 |
| 60075783 | 36 |
| 60076302 | 36 |
| 60077361 | 10 |
| 60077478 | 20 |
| 60078314 | 34, 36 |
| 60078383 | 36 |
| 60078475 | 32 |
| 60078507 | 34 |

| Part No. | Page(s) |
|----------|---------|
| 60078508 | 34 |
| 60078510 | 38 |
| 60078516 | 34 |
| 60078517 | 34 |
| 60078519 | 34 |
| 60078530 | 34 |
| 60078531 | 24 |
| 60078537 | 34 |
| 60078539 | 36 |
| 60078940 | 38 |
| 60079419 | 28, 32 |
| 60079472 | 36 |
| 60079612 | 8 |
| 60080240 | 8 |
| 60080871 | 24 |
| 60081489 | 46 |
| 60081519 | 46 |
| 60083598 | 32 |
| 60084455 | 22 |
| 60084479 | 36 |
| 60084840 | 22 |
| 60086446 | 42 |
| 60086467 | 36 |
| 60086575 | 32 |
| 60087522 | 26 |

| Part No. | Page(s) |
|----------|---------|
| 60088416 | 42 |
| 60090139 | 28, 32 |
| 60090207 | 36 |
| 60090830 | 34 |
| 60091424 | 42 |
| 60091427 | 42 |
| 60091949 | 36 |
| 60091950 | 32 |
| 60092440 | 28, 32 |
| 60092630 | 22 |
| 60092766 | 50 |
| 60093006 | 10 |
| 60093007 | 20 |
| 60093008 | 22 |
| 60093009 | 28, 32 |
| 60093452 | 46 |
| 60093500 | 26 |

2. Introduction

PROPRIETARY
OF FOSS

| | |
|---|---|
| X On Site Maintenance | |
| <input type="checkbox"/> Remote Maintenance | |
| | |
| PM Interval | 12 Months |
| PM Kit (P/N) | 6010 9491 |
| Valid From | 2026-01-15 |
| Additional P/N need | 1 × 6009 3021 (FTIR Performance Standard 753) |
| | |
| Instrument Serial No.: | |
| No. of Samples Analyzed: | |
| Date: | |

A maintenance protocol provides a systematic and effective approach to ensure the reliable operation of a specific instrument type. The recommended Preventive Maintenance (PM) interval depends on operational conditions and is guided by our extensive experience in manufacturing and servicing analytical instruments.

In addition to sample throughput, environmental factors must also be taken into account. Challenging conditions - such as elevated ambient temperatures, high humidity, or excessive dust - can significantly reduce component lifetimes and shorten both maintenance and replacement intervals.

The maintenance frequency of the MilkoScan™ FT3 is provided in the table above.

2.1. Maintenance Procedures

1. Pre-verification of Instrument before PM replacement

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.3):



| No. | Instrument Performance Test | OK |
|----------------------|--|--------------------------|
| 1 | Perform the Instrument Verification before the PM starts | <input type="checkbox"/> |
| 2 | Perform a Clean and a Zero Setting | <input type="checkbox"/> |
| Time required | | 0.2 hour |

Table 1 - Pre-verification

2. Removal and Replacement

See the Spare Parts Manual (SP 6008 5771) and procedure in the Service Manual (SM 6008 5772 chapter 3.3.4):

| Part / Module | Page / Pos | Action | Qty | P/N | OK |
|---|------------|--------------------------|-----|--------------|--------------------------|
| Membrane, ø12 for Back Pressure Valve | 23/6 | Replace | 1 | 6005 3831 | <input type="checkbox"/> |
| Cuvette MilkoScan FT3 SP (in plastic box) | 20/3 5 | Replace/ Hand over | 1 | 6010 0192 | <input type="checkbox"/> |
| O-ring Ø003.00xØ1.00 for Inline Filter | 20/2 | Replace | 1 | 6001 5107 | <input type="checkbox"/> |
| Inline Filter, Mesh Filter | 20 | Check/ve rify | 1 | | <input type="checkbox"/> |
| O-ring ø002.50/ø1.00 FPM 80 | 30/18 | Replace | 1 | 6001 5423 | <input type="checkbox"/> |
| O-ring ø003.70xø1.20 FPM 70 | 30/17 | Replace | 1 | 6010 9865 | <input type="checkbox"/> |
| Tube ø1/8"/ø01.60 PTFE L=0130mm assy | 40/5 6 | Replace | 1 | 6007 0955 | <input type="checkbox"/> |
| Tube ø01.59/ø03.18 PTFE L=0725 assy | 40/5 3 | Replace | 1 | 6007 4464 | <input type="checkbox"/> |
| Tube Ø01.59/Ø03.18 PTFE L=0415 assy | 40/5 4 | Replace | 1 | 6007 5536 | <input type="checkbox"/> |
| Tube ø01.59/ø03,18 PTFE L=0145 assy 1 | 40/5 50 | Replace | 1 | 6009 1427 | <input type="checkbox"/> |
| Tube øø 008/016 Silicone as Waste tube | 6/4 | Replace | 1 | 6007 8292 | <input type="checkbox"/> |
| Repair Kit Syringe Piston Tip | 30/B | Replace | 1 | 6010 8261 | <input type="checkbox"/> |
| Motor step 43×43mm, Screw L=2 in cpl 2 | 30/12 | Grease | 2 | 6009 0139 | <input type="checkbox"/> |

| Part / Module | Page / Pos | Action | Qty | P/N | OK |
|---|------------|---------|-----|--------------|--------------------------|
| Repair Kit Valve Flow house ø0.8, FKM | 30/2 2 | Replace | 2 | 6010 9532 | <input type="checkbox"/> |
| Repair Kit Syringe Piston Tip | 26/B | Replace | 1 | 6010 8261 | <input type="checkbox"/> |
| Motor step 43×43mm, Screw L=2 in cpl 2 | 26/11 | Grease | 2 | 6009 0139 | <input type="checkbox"/> |
| Repair Kit Valve Flow house ø1.2, FKM | 26/B | Replace | 1 | 6010 9531 | <input type="checkbox"/> |
| Label sheet, Clean and Zero for Container Lids and Containers | | Replace | 1 | 6008 8049 | <input type="checkbox"/> |
| Time required | | | | | 3 hours |

Table 2 - Remove/Replace

1. Only if instrument is equipped with Conductivity Module
2. Motor Spindle should be lightly greased with Kryton 206 (P/N 6007 9634 Lubricant 57 g).

3. Post-verification of Instrument after PM

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.5):

| No. | Instrument Performance Test | Batch number | OK |
|----------------------|--|--------------|--------------------------|
| 1 | Perform the Instrument Verification by the end of PM | | <input type="checkbox"/> |
| 2 | Perform a Clean and a Zero Setting | N/A | <input type="checkbox"/> |
| Time required | | | 0.2 hour |

Table 3 - Post-verification

4. Instrument Check

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.6):

| No | Action | OK |
|----------------------|--|--------------------------|
| . | | |
| 1 | Check Diagnostic center for any detected events (no events accepted) | <input type="checkbox"/> |
| 2 | Check (view log) and create an "event log report" in Instrument performance reports | <input type="checkbox"/> |
| 3 | Perform "Back up the entire database" | <input type="checkbox"/> |
| 4 | Perform "Clean up data" | <input type="checkbox"/> |
| 5 | Redo "Back up the entire database" | <input type="checkbox"/> |
| 6 | Check "Flow system check" (no warnings/errors accepted) | <input type="checkbox"/> |
| 7 | Check "Sample pipette leak check" (no warnings/errors accepted) | <input type="checkbox"/> |
| 8 | Check "IR system check - measure with no cuvette flow" (no warnings/errors accepted) | <input type="checkbox"/> |
| 9 | Check "Flow system leak diagnostic tool" (no warnings/errors accepted) | <input type="checkbox"/> |
| 10 | Export Wear log at "Save service data to file". Send to Global Help Desk | <input type="checkbox"/> |
| Time required | | 0.5 hour |

Table 4 - Instrument Check

5. Software

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.7):

| Action | Version Number | OK |
|--|----------------|--------------------------|
| Currently installed version of Nova | | <input type="checkbox"/> |
| Perform an upgrade to the latest version (if needed) | | <input type="checkbox"/> |
| Time required | | 0.5 hour |

Table 5 - Software verification

6. Instrument Performance Verification

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.8):

| No. | Instrument Performance Test | OK |
|----------------------|--|--------------------------|
| 1 | Perform a Zero Setting | <input type="checkbox"/> |
| 2 | Perform a Clean and a Zero Setting | <input type="checkbox"/> |
| 3 | Repeatability test on homogenized whole milk | <input type="checkbox"/> |
| 4 | Carry over test | <input type="checkbox"/> |
| 5 | Volume check | <input type="checkbox"/> |
| 6 | Perform a Clean and a Zero | <input type="checkbox"/> |
| Time required | | 0.5 hour |

Table 6 - Instrument verification

| | |
|--|--|
| <p>Note: The total time is an approximate time. Depending on e.g. the current status and the accessibility of the instrument, this time might vary. Please round off to nearest quarter of an hour.</p> | <p>Total Time [Hours] 5</p> |
|--|--|

| | |
|--------------------------|------------|
| Customer representative: | Signature: |
| FOSS representative: | Signature: |

2.1.1. 1. Pre-verification of Instrument before PM replacement

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.3):

| No. | Instrument Performance Test | OK |
|----------------------|--|--------------------------|
| 1 | Perform the Instrument Verification before the PM starts | <input type="checkbox"/> |
| 2 | Perform a Clean and a Zero Setting | <input type="checkbox"/> |
| Time required | | 0.2 hour |

Table 1 - Pre-verification

2.1.2. 2. Removal and Replacement

See the Spare Parts Manual (SP 6008 5771) and procedure in the Service Manual (SM 6008 5772 chapter 3.3.4):

| Part / Module | Page / Pos | Action | Qty | P/N | OK |
|---|------------|--------------------------|-----|--------------|-----|
| Membrane, ø12 for Back Pressure Valve | 23/6 | Replace | 1 | 6005 3831 | [] |
| Cuvette MilkoScan FT3 SP (in plastic box) | 20/3 5 | Replace/ Hand over | 1 | 6010 0192 | [] |
| O-ring Ø003.00xØ1.00 for Inline Filter | 20/2 | Replace | 1 | 6001 5107 | [] |
| Inline Filter, Mesh Filter | 20 | Check/ve rify | 1 | | [] |
| O-ring ø002.50/ø1.00 FPM 80 | 30/18 | Replace | 1 | 6001 5423 | [] |
| O-ring ø003.70xø1.20 FPM 70 | 30/17 | Replace | 1 | 6010 9865 | [] |
| Tube ø1/8"/ø01.60 PTFE L=0130mm assy | 40/5 6 | Replace | 1 | 6007 0955 | [] |
| Tube ø01.59/ø03.18 PTFE L=0725 assy | 40/5 3 | Replace | 1 | 6007 4464 | [] |
| Tube Ø01.59/Ø03.18 PTFE L=0415 assy | 40/5 4 | Replace | 1 | 6007 5536 | [] |
| Tube ø01.59/ø03,18 PTFE L=0145 assy 1 | 40/5 50 | Replace | 1 | 6009 1427 | [] |
| Tube øø 008/016 Silicone as Waste tube | 6/4 | Replace | 1 | 6007 8292 | [] |
| Repair Kit Syringe Piston Tip | 30/B | Replace | 1 | 6010 8261 | [] |
| Motor step 43×43mm, Screw L=2 in cpl 2 | 30/12 | Grease | 2 | 6009 0139 | [] |

| Part / Module | Page / Pos | Action | Qty | P/N | OK |
|---|------------|---------|-----|--------------|------------|
| Repair Kit Valve Flow house ø0.8, FKM | 30/2 2 | Replace | 2 | 6010 9532 | [] |
| Repair Kit Syringe Piston Tip | 26/B | Replace | 1 | 6010 8261 | [] |
| Motor step 43×43mm, Screw L=2 in cpl 2 | 26/11 | Grease | 2 | 6009 0139 | [] |
| Repair Kit Valve Flow house ø1.2, FKM | 26/B | Replace | 1 | 6010 9531 | [] |
| Label sheet, Clean and Zero for Container Lids and Containers | | Replace | 1 | 6008 8049 | [] |
| Time required | | | | | 3 hours |

Remove/Replace

1 - Only if instrument is equipped with Conductivity Module.

2 - Motor Spindle should be lightly greased with Kryton 206 (P/N 6007 9634 Lubricant 57 g).

2.1.3. 3. Post-verification of Instrument after PM

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.5):

| No. | Instrument Performance Test | Batch number | OK |
|----------------------|--|--------------|--------------------------|
| 1 | Perform the Instrument Verification by the end of PM | | <input type="checkbox"/> |
| 2 | Perform a Clean and a Zero Setting | N/A | <input type="checkbox"/> |
| Time required | | | 0.2 hour |

Table 3 - Post-verification

2.1.4. 4. Instrument Check

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.6):

| No | Action | OK |
|----------------------|--|--------------------------|
| . | | |
| 1 | Check Diagnostic center for any detected events (no events accepted) | <input type="checkbox"/> |
| 2 | Check (view log) and create an "event log report" in Instrument performance reports | <input type="checkbox"/> |
| 3 | Perform "Back up the entire database" | <input type="checkbox"/> |
| 4 | Perform "Clean up data" | <input type="checkbox"/> |
| 5 | Redo "Back up the entire database" | <input type="checkbox"/> |
| 6 | Check "Flow system check" (no warnings/errors accepted) | <input type="checkbox"/> |
| 7 | Check "Sample pipette leak check" (no warnings/errors accepted) | <input type="checkbox"/> |
| 8 | Check "IR system check - measure with no cuvette flow" (no warnings/errors accepted) | <input type="checkbox"/> |
| 9 | Check "Flow system leak diagnostic tool" (no warnings/errors accepted) | <input type="checkbox"/> |
| 10 | Export Wear log at "Save service data to file". Send to Global Help Desk | <input type="checkbox"/> |
| Time required | | 0.5 hour |

Table 4 - Instrument Check

2.1.5. 5. Software

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.7):

| Action | Version Number | OK |
|--|----------------|--------------------------|
| Currently installed version of Nova | | <input type="checkbox"/> |
| Perform an upgrade to the latest version (if needed) | | <input type="checkbox"/> |
| Time required | | 0.5 hour |

Table 5 - *Software verification*

2.1.6. 6. Instrument Performance Verification

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.8):

| No. | Instrument Performance Test | OK |
|----------------------|--|--------------------------|
| 1 | Perform a Zero Setting | <input type="checkbox"/> |
| 2 | Perform a Clean and a Zero Setting | <input type="checkbox"/> |
| 3 | Repeatability test on homogenized whole milk | <input type="checkbox"/> |
| 4 | Carry over test | <input type="checkbox"/> |
| 5 | Volume check | <input type="checkbox"/> |
| 6 | Perform a Clean and a Zero | <input type="checkbox"/> |
| Time required | | 0.5 hour |

Table 6 - Instrument verification

Note: The total time is an approximate time. Depending on e.g. the current status and the accessibility of the instrument, this time might vary. Please round off to nearest quarter of an hour.

**Total
Time
[Hours]**
5

3. **Maintenance Procedures (w/Index)**

3.1. 1. Pre-verification of Instrument before PM replacement

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.3):

| No. | Instrument Performance Test | OK |
|----------------------|--|--------------------------|
| 1 | Perform the Instrument Verification before the PM starts | <input type="checkbox"/> |
| 2 | Perform a Clean and a Zero Setting | <input type="checkbox"/> |
| Time required | | 0.2 hour |

Table 1 - Pre-verification

3.2. 2. Removal and Replacement

See the Spare Parts Manual (SP 6008 5771) and procedure in the Service Manual (SM 6008 5772 chapter 3.3.4):

| Part / Module | Page / Pos | Action | Qty | P/N | OK |
|---|------------|--------------------------|-----|--------------|-----|
| Membrane, ø12 for Back Pressure Valve | 23/6 | Replace | 1 | 6005 3831 | [] |
| Cuvette MilkoScan FT3 SP (in plastic box) | 20/3 5 | Replace/ Hand over | 1 | 6010 0192 | [] |
| O-ring Ø003.00xØ1.00 for Inline Filter | 20/2 | Replace | 1 | 6001 5107 | [] |
| Inline Filter, Mesh Filter | 20 | Check/ve rify | 1 | | [] |
| O-ring ø002.50/ø1.00 FPM 80 | 30/18 | Replace | 1 | 6001 5423 | [] |
| O-ring ø003.70xø1.20 FPM 70 | 30/17 | Replace | 1 | 6010 9865 | [] |
| Tube ø1/8"/ø01.60 PTFE L=0130mm assy | 40/5 6 | Replace | 1 | 6007 0955 | [] |
| Tube ø01.59/ø03.18 PTFE L=0725 assy | 40/5 3 | Replace | 1 | 6007 4464 | [] |
| Tube Ø01.59/Ø03.18 PTFE L=0415 assy | 40/5 4 | Replace | 1 | 6007 5536 | [] |
| Tube ø01.59/ø03,18 PTFE L=0145 assy 1 | 40/5 50 | Replace | 1 | 6009 1427 | [] |
| Tube øø 008/016 Silicone as Waste tube | 6/4 | Replace | 1 | 6007 8292 | [] |
| Repair Kit Syringe Piston Tip | 30/B | Replace | 1 | 6010 8261 | [] |
| Motor step 43×43mm, Screw L=2 in cpl 2 | 30/12 | Grease | 2 | 6009 0139 | [] |

| Part / Module | Page / Pos | Action | Qty | P/N | OK |
|---|------------|---------|-----|--------------|------------|
| Repair Kit Valve Flow house ø0.8, FKM | 30/2 2 | Replace | 2 | 6010 9532 | [] |
| Repair Kit Syringe Piston Tip | 26/B | Replace | 1 | 6010 8261 | [] |
| Motor step 43×43mm, Screw L=2 in cpl 2 | 26/11 | Grease | 2 | 6009 0139 | [] |
| Repair Kit Valve Flow house ø1.2, FKM | 26/B | Replace | 1 | 6010 9531 | [] |
| Label sheet, Clean and Zero for Container Lids and Containers | | Replace | 1 | 6008 8049 | [] |
| Time required | | | | | 3 hours |

Remove/Replace

1 - Only if instrument is equipped with Conductivity Module.

2 - Motor Spindle should be lightly greased with Kryton 206 (P/N 6007 9634 Lubricant 57 g).

3.3. 3. Post-verification of Instrument after PM

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.5):

| No. | Instrument Performance Test | Batch number | OK |
|----------------------|--|--------------|--------------------------|
| 1 | Perform the Instrument Verification by the end of PM | | <input type="checkbox"/> |
| 2 | Perform a Clean and a Zero Setting | N/A | <input type="checkbox"/> |
| Time required | | | 0.2 hour |

Table 3 - Post-verification

3.4. 4. Instrument Check

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.6):

| No | Action | OK |
|----------------------|--|--------------------------|
| . | | |
| 1 | Check Diagnostic center for any detected events (no events accepted) | <input type="checkbox"/> |
| 2 | Check (view log) and create an "event log report" in Instrument performance reports | <input type="checkbox"/> |
| 3 | Perform "Back up the entire database" | <input type="checkbox"/> |
| 4 | Perform "Clean up data" | <input type="checkbox"/> |
| 5 | Redo "Back up the entire database" | <input type="checkbox"/> |
| 6 | Check "Flow system check" (no warnings/errors accepted) | <input type="checkbox"/> |
| 7 | Check "Sample pipette leak check" (no warnings/errors accepted) | <input type="checkbox"/> |
| 8 | Check "IR system check - measure with no cuvette flow" (no warnings/errors accepted) | <input type="checkbox"/> |
| 9 | Check "Flow system leak diagnostic tool" (no warnings/errors accepted) | <input type="checkbox"/> |
| 10 | Export Wear log at "Save service data to file". Send to Global Help Desk | <input type="checkbox"/> |
| Time required | | 0.5 hour |

Table 4 - Instrument Check

3.5. 5. Software

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.7):

| Action | Version Number | OK |
|--|----------------|--------------------------|
| Currently installed version of Nova | | <input type="checkbox"/> |
| Perform an upgrade to the latest version (if needed) | | <input type="checkbox"/> |
| Time required | | 0.5 hour |

Table 5 - *Software verification*

3.6. 6. Instrument Performance Verification

See the procedure in the Service Manual (SM 6008 5772 chapter 3.3.8):

| No. | Instrument Performance Test | OK |
|----------------------|--|--------------------------|
| 1 | Perform a Zero Setting | <input type="checkbox"/> |
| 2 | Perform a Clean and a Zero Setting | <input type="checkbox"/> |
| 3 | Repeatability test on homogenized whole milk | <input type="checkbox"/> |
| 4 | Carry over test | <input type="checkbox"/> |
| 5 | Volume check | <input type="checkbox"/> |
| 6 | Perform a Clean and a Zero | <input type="checkbox"/> |
| Time required | | 0.5 hour |

Table 6 - Instrument verification

Note: The total time is an approximate time. Depending on e.g. the current status and the accessibility of the instrument, this time might vary. Please round off to nearest quarter of an hour.

**Total
Time
[Hours]**
5

4. Untitled
