

Dental Accessory Preparation Checklist

INTRODUCTION

This document lists the actions that must be taken to unbox and organize your accessories before familiarization training. Please make sure all the steps below are completed before your scheduled training.

1. Unwrap blue mats and tape down to tables/workspace with painter's tape.
2. Unbox IPA and acetone bottles.
 - a. Fill each bottle up to marked line with the appropriate solvent.
 - b. Place remainder of the solvent in the flammables cabinet for quick access.
3. Unbox DPR 10 resin and place in storage cabinet.
4. Unwrap orbital shaker and connect to power.
 - a. Cut a blue mat to size and add to top of orbital shaker to absorb spillage.
5. Ensure the following items are present:

a. Two plastic trays



b. Nitrile gloves



c. Blue towels



d. Aluminum foil



e. Dead blow hammer



f. Yellow scraper



6. Cover both trays with a layer of foil.
7. Unbox the cassette pour stand and set on bench.
8. Place the following items by the cassette pour stand:

a. Funnels



b. Paint strainers



c. Spatula



d. Blue towels

